


Notch rejection filters

TopNotch™


Chroma Technology manufactures high-performance optical interference filters using sputtering and ion beam technologies. Covering a spectral range from 200nm- 3000nm, our wide-ranging capabilities enable us to offer a broad array of filter products offering superior durability and longevity, and spanning many industrial applications.




Best-in-Class Performance

- Extra-wide, high-transmission range:
350 - 1600nm
- Narrow notch:
approx. 3% of notch CWL
- High blocking level:
>OD6
- Wide blocking range:
>6nm wavelength range

Narrow band, square rejection



Wider range of high, OD6 blocking



ZET488 TopNotch

ZET532 TopNotch

ZET594 TopNotch

ZET642 TopNotch

ZET785 TopNotch

CHROMA TECHNOLOGY CORP®

Contact: sales@chroma.com or +1-800-824-7662 • +1-802-428-2500 • chroma.com

